

Cat. No. C15-6052-02 Vacuum buzzer

SPECIFICATIONS

Vacuum buzzer: (Vibration observation type)

The buzzer is driven by a single size AA battery. You can see a small styrene foam ball in contact with the diaphragm of the buzzer jumps around and the buzzer sounds. (Battery is not included.)

Size: Approx. $30 \times 34 \times 74$ mm

Weight: Approx. 29 g (Battery is not included.)

Experiment using the vacuum buzzer

(See Figs. 2 and 3.)

- (1) Attach the small rubber band (one) included in the product into the back of the vacuum cover with adhesive tape as shown in Fig. 2.



Fig. 2

Rubber band
Neo sponge
packing
Adhesive tape
Vacuum cover

- (2) Then, connect a size D battery to the vacuum buzzer. The buzzer sounds and the small styrene foam ball in the acrylic pipe jumps around as the buzzer vibrates. (The size D battery is not included.)
Note: Turn off the vacuum buzzer simply by putting paper between the + pole of the battery and connector and turn on the buzzer by removing it.
- (3) Hang the sounding buzzer on the rubber band attached to the vacuum cover and put them into the vacuum container carefully to avoid the buzzer coming into contact with the wall of the container.
- (4) Make sure that the hung buzzer generates a quiet sound. The styrene foam ball jumps around and the buzzer sounds.





Fig. 3

The styrene foam

- (5) Discharge air from the container. In this condition, the ball of styrene foam vibrates and you cannot hear the buzzer sound.
- (6) Tilt the vacuum container. Visually, you can see that the buzzer keeps sounding.



In this condition, the ball of styrene foam vibrates and you cannot hear the buzzer sound.

Note: Sound is transmitted in the air, water, and solid matters but is not transmitted in a vacuum. In this experiment, the container does not reach a complete vacuum so you can hear a quiet sound in a silent environment.

- (7) Loosen the pinch cock. Air is let into the container and the buzzer sound is heard again.

NARIKA CORPORATION

<http://www.rika.com>

Head office Soto-Kanda 5-3-10, Chiyoda-ku, Tokyo, 101-0021, Japan
TEL +81-3-3833-0741 (main) FAX +81-3-3836-1725

Osaka office Nagata Nakatsu Bld. 5th floor, Oyodo Naka 1-4-16, Kita-ku, Osaka-shi, Osaka, 531-0076, Japan
TEL +81-6-6451-3986 FAX +81-6-6451-3925

Fukuoka office Higashi Hakata Ace Bld. 7th floor, Hie-cho 2-7, Hakata-ku, Fukuoka-shi, Fukuoka, 812-0014, Japan
TEL +81-92-432-6888 FAX +81-92-432-7388

Inquiry about the products ...

Science Tool Counselor Office 0120-700-746